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The use
of diagnostic
systems
in professional
sports
and fitness
training


Timur Bestavishvili

Sport has become a lifestyle choice for millions of people and one of the leading industries. Like other rapidly developing industries, that of sport is in urgent need of instrumentation. This term refers not only to sports facilities and equipment or result-oriented techniques. A major role is increasingly played by means of monitoring the health status during sports activities. In this area the importance of rapid diagnosis and correction equipment such as the Omega M System can hardly be overstated

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The modern professional sports coaching techniques are so highly developed that they make it possible to train virtually 24 hours a day with continuously improving performance results. There is only one problem: that of restoring the body to the original state of physical fitness after a routine training session. The problem is therefore not one of sustaining exercise but one of recuperating after it.

*A number of effective recuperation techniques have been developed, which fall into four groups:
– physiological techniques based on physiotherapy such as sauna treatment and various types of massage;*

*– methodological techniques based on alternative types of physical activity stimulating recuperation processes;
– pharmacological techniques based on the administration of a group of compensatory drugs (it should be pointed out that a number of such drugs can be plausibly classified as doping agents given their adverse long-term effects);
– psychophysiological techniques based on control over the person's mental state (creating a certain state of mind targeting recuperation and relaxation).*

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Obviously, diagnosis should also be given a thought to. Rapid diagnostic test results are certainly insufficient to form a definitive diagnosis, but such a diagnosis is not needed in this case. Should the sports therapist detect even a suggestion of something being wrong with the athlete's body, he or she can always refer the athlete to the appropriate specialist for advanced diagnostic testing. That is what makes it necessary to employ diagnostic equipment that would facilitate this kind of work.

Such diagnostic equipment should have the following properties:

- be easy to use;*
- enable rapid (minutes-long) basic physical examination;*
- have an informative interface, easily usable by physicians, coaches, and athletes, including amateurs;*
- support the processing of all relevant medical data and their transmission to any specialist, both electronically and in printed form;*
- enable follow-up supervision, making it possible to detect even the slightest changes in the athlete's condition;*
- enable the overall assessment of the athlete's form and fitness;*

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– support the correction of the athlete's condition by a specialist on the basis of the data obtained.

It is especially important to take rapid diagnostic results into account when prescribing a course of treatment. A rapid analysis system is certainly unable to provide a proper diagnosis; however, it can clearly indicate the existence of a problem. After this a goal-oriented search for the origin of the disorder can be conducted using more sophisticated and accurate diagnostic methods, whereupon the treatment can be continued or modified if necessary. The methodological bases of rapid diagnosis include not only a mechanism for


implementing the method, but also a task-oriented database for recording and tracking examination results. The database supports medical record-keeping both in the electronic filing format, familiar to any physician, and a less familiar format, which is, however, more convenient to use when dealing with specific tasks. Naturally, the electronic mode of operation does not exclude the hard copy option. The equipment for comprehensive computer-aided functional diagnosis

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of the human body includes the Omega M and Omega C Systems as well as the auxiliary Neuron and Crab Systems. Their interfaces are easy to use and amateur-friendly. The numerous explanatory notes accessible through the help menu, as well as the convenient electronic record-keeping (filing) system, will safely guide you through the options available. The Electronic Patient Medical Record feature makes it possible to record and store test results, information about the treatments prescribed, specialists' diagnoses, and other data normally stored in a 'paper' filing system. The Omega System uses electromagnetic pulses to measure

certain functional parameters of the body and specifically to assess deviations from the physiological optimum. Moreover, the system provides a basis for conclusions concerning the nature of the condition in question – whether it is a heart, nervous, or endocrine disorder, etc. All these data can be obtained by I-lead ECG analysis (using a pair of wrist electrodes) in 5 to 6 minutes. An advantage of Omega software is that it uses analysis results to calculate a corrective forced respiration treatment algorithm. Cardiac function is well known to be influenced to a certain extent by blood carbon dioxide levels. An inspiration-expiration-pause

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
rhythm calculated by Omega and displayed on the computer screen helps considerably to improve the patient's condition in about four minutes. The effect of such treatment is likely to be transient, though. Another Omega-based strategy is available to achieve a lasting effect. The diagnostic system makes it possible to administer intensive and effective corrective treatment based on balancing the multiple and complex electromagnetic rhythms of the body. Such treatment employs a personalized audiovisual correction program created using the patient's diagnostic test results and written to a DVD. The method has proved highly effective, especially as a means of weight control and exercise

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
management. Reports on this are available in the literature on the subject.

Thus, the high-performance real-time diagnostic system under consideration enables corrective treatment and meets all the above requirements. Omega M software is intended for use by physicians, and especially sports therapists. Omega C offers virtually the same functionality, but its interface is adapted for use by coaches and athletes. Neuron is a full-scale 12-lead

digital ECG monitor, and Crab is a system for monitoring seven or more patients simultaneously. The entire set of equipment as well as its operation principles is patented and certified by numerous regulatory authorities. Given the current trends in the sports services market, a medical facility equipped with a diagnostic system is an integral image-building part of an up-to-date sports club or spa salon regardless of the format it operates in. ■

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© 2009 Jupiterimages Corporation